



Accuflow[®] Thala Pump



vTitan Corporation Private Limited

New Woodlands Building, Estancia IT Park,
Plot No. 140,151 Vallancheri Village,
Chengalpattu District,
Tamil Nadu - 603202, India.
Toll Free Number: 18005721664

www.vtitan.com

sales@vtitan.com

©2026 vTitan Corporation Pvt. Ltd.





980-1062 Rev 10

A Helping Hand to Fight Thalassemia...


Why Thala Pump :

Thalassemia, an inherited blood disorder characterized by less haemoglobin often requires frequent blood transfusions, resulting in iron overload. Thala pump infuses chelating agents to remove the excess iron.

Product Features:

-  Compact and Easy to use
-  Simple priming procedure
-  Lock functionality
-  Visual and audio alarms

3-Step Operation



Product Specifications:

Dimension (in mm)	172 x 36 x 51 (L x W x H)
Weight (in grams)	185
Compatible syringe sizes	5, 10, 20 mL
Display	0.96 inch OLED
Flow rate range	0.1 to 70.0 mL/hr
Accuracy	± 2%
Priming function	Available
Battery type	9V Li-Ion battery (Rechargeable)
Battery life	40 hrs
Alarms	Occlusion, Syringe near empty, Syringe empty, Syringe disengaged, Battery low, Battery depleted, Unit malfunction
Accessories	Pouch, charger, straps
Warranty	1 year